


LED OPTOTYPE OTTOTIPI LED



BFSL 27.67 FL

BFSL 27.67 EL

BFSL 27.67 DL

BFSL 27.67 CL

BFSL 27.67 BL

BFSL 27.67 AL


CONFIGURATION

- Chassis: body composed of anti scratch anodized aluminium sections
- Light Source: LED technology, the most advanced technology applied to lighting systems with 5000 cd/mq. Temperature of colour 7200° K
- Diffusing Panel: white acrylic sheet of 3 mm thickness
- Switch: ON/OFF
- Power supply: 230V; 50 Hz
- Compliant with essential requirements of the EC 2007/47 as Medical Device Class I

CONFIGURAZIONE

- Chassis: realizzato in profilato di alluminio anodizzato antigraffio
- Fonte Luminosa: impiego della tecnologia LED, la più avanzata tecnologia applicata ai sistemi di illuminazione con 5000 cd/mq. Temperatura colore 7200° K
- Pannello Diffusore: lastra acrilica bianco opalino, spessore 3 mm
- Interruttore: ON/OFF
- Alimentazione: 230V; 50 Hz
- Conforme ai requisiti essenziali della EC 2007/47 in qualità di Dispositivo Medico di Classe I

Code/Codice	Type/Modello	Dimensions/Dimensioni W x D x H [cm]	Weight/Peso [kg]
BFSL 27,67 AL	Mixed decimal <i>Decimale mista</i>	3 mt	34 x 3 x 72
BFSL 27,67 BL	Monoyer <i>Monoyer</i>	3 mt	35 x 3 x 72
BFSL 27,67 CL	Children <i>Bambini</i>	3 mt	36 x 3 x 72
BFSL 27,67 DL	Mixed decimal <i>Decimale mista</i>	5 mt	37 x 3 x 72
BFSL 27,67 EL	Armagnac <i>Armagnac</i>	5 mt	38 x 3 x 72
BFSL 27,67 FL	Children <i>Bambini</i>	5 mt	39 x 3 x 72



weiko[®]
health base

WEIKO - BIEFFE ITALIA s.r.l.

Sede legale:
Viale Lincoln, 178
81100 Caserta
info@weiko.com
www.weiko.com

Stabilimento:
Consorzio Sinè
Zona ASI Aversa Nord
81032 Carinaro (CE)
Tel. +39 081 8654662
Fax: +39 081 8654723

www.weiko.com