

Weiko®

www.weiko.com

Technical Innovations and Design



Optotypes

“Luminous” Optotypes

- **Chassis:** body composed of oven-fired epoxy powder coated aluminium sections with nylon junctions
- **Light Source:** straight daylight fluorescent tubes 18W
- **Diffusing Panel:** acrylic sheet of 3 mm thickness
- **Power Supply:** 230V; 50/60 Hz
- **Switch:** ON/OFF bipolar switch with two

- **pole fused plug.**
- **Cable:** electric cable made of pre-insulated terminals in self extinguishing PVC
- **Components:** all components UL - DVE - IMQ approved. All optotypes correspond to EC regulations.
- **Conformity:** the appliance complies with the essential requirements of the

Low Voltage Directive 73/23 EEC (CEI EN 60598-1 and upgradings A1/A12/A13; CEI EN 60598-2-4; CEI EN 60601-1 and upgradings A1/A2/A11/A12/A13) and with the essential requirements of the Electro Magnetic Compatibility 89/336 EEC (CEI EN 55015/97 and upgradings A1/A2; CEI EN 61547; CEI EN 61000-3-2/95 and upgradings A1/A2).

Luminous Types	Code	N° switches	Viewing area (cm)	External dimensions (cm)			Net weight (Kg)	Temperature of colour	Luminous flux (lm)	Chromatic Index	Total Power (Watt)
				Length	Height	Depth					
Mixed decimal 3 mt	BF AL 27.67AL	1	24x62	29,5	67,5	12	6	6500	700	63	18
Monoyer 3 mt	BF AL 27.67BL	1	24x62	29,5	67,5	12	6	6500	700	63	18
Children 3 mt	BF AL 27.67CL	1	24x62	29,5	67,5	12	6	6500	700	63	18
Mixed decimal 5 mt	BF AL 27.67DL	1	24x62	29,5	67,5	12	6	6500	700	63	18
Armagnac 5 mt	BF AL 27.67EL	1	24x62	29,5	67,5	12	6	6500	700	63	18
Children 5 mt	BF AL 27.67FL	1	24x62	29,5	67,5	12	6	6500	700	63	18
		BF-RC	Remote control: available for all luminous optotypes - switch ON/OFF powerful up to 30 mt								

“Non Luminous” Optotypes

- **Chassis:** body composed of oven-fired epoxy powder coated aluminium sections with nylon junctions

- **Panel:** MDF sheet of 3 mm thickness

Non Luminous Types	Code	Viewing area (cm)	External dimensions (cm)			Net weight (Kg)
			Length	Height	Depth	
Mixed decimal 3 mt	BF AL 27.67A	24x61	29	66	2,5	2
Monoyer 3 mt	BF AL 27.67B	24x61	29	66	2,5	2
Children 3 mt	BF AL 27.67C	24x61	29	66	2,5	2
Mixed decimal 5 mt	BF AL 27.67D	24x61	29	66	2,5	2
Armagnac 5 mt	BF AL 27.67E	24x61	29	66	2,5	2
Children 5 mt	BF AL 27.67F	24x61	29	66	2,5	2

Acronym Legend

Code	Standard models
A / AL	Mixed decimal 3 mt
B / BL	Monoyer 3 mt
C / CL	Children 3 mt
D / DL	Mixed decimal 5 mt
E / EL	Armagnac 5 mt
F / FL	Children 5 mt
Code	On request
G / GL	Landolt 3 mt
H / HL	Landolt 5 mt
I / IL	Wecker 5 mt
L / LL	Intern. "ESS" 5 mt

Weiko Italia

Viale Lincoln, 152
81100 - Caserta
info@weiko.com
www.weiko.com

Via Palma, 284
80040 - Striano (NA)
Tel. +39 081 649857
Fax +39 081 643324

